

## Disinfectant cleaning solution



Powerfully effective, independently tested to British and European standards, and suitable for use on multiple contact surfaces.

### Performance focus

Cutting-edge research and development, combined with independent testing, ensures **Bioguard disinfectant cleaning solution** is truly effective for use across the whole spectrum of professional healthcare environments.

Our intensive focus on performance delivers a product which is:

- POWERFULLY EFFECTIVE** - Independently tested to British and European Standards BSEN1276. Effective against 1000's of viruses, bacteria, mycobacterium and spores, including MRSA, Tuberculosis, Hepatitis C, Salmonella, Avian Flu (H5N1) and Swine Flu (H1N1). For more details, please refer to our Independent Test and Safety Information CD.
- NON-HAZARDOUS** - Conforming to the latest health and safety legislation. For more details, please refer to our Independent Test and Safety Information CD.
- USABLE ON VIRTUALLY ALL CONTACT SURFACES** - Including plastic, stainless steel, bedding and uniforms, without damaging or degrading them; one quick, hassle-free application is all that's required.
- SUPERIOR IN IT'S CLEANSING ACTION** - independently proven to out-perform other leading European products.
- IDEAL FOR HIGH-RISK ENVIRONMENTS** - The professionals' choice since the early 1990s; widely used in clinical facilities, including mortuaries and testing labs.
- HIGHLY ECONOMIC** - Remains effective at higher dilution rates, reducing the amount of applied product solution used.
- DUAL ACTION** - Combined cleaning and disinfecting action eliminates the need for separate products and procedures.
- ENVIRONMENTALLY FRIENDLY** - Proven by independent audits to be non-toxic, biodegradable and non-accumulative in the eco-system. No special procedures are required for product disposal or in case of spillage.
- USER FRIENDLY** - Formulated to prevent skin irritation or unusual sensitivity.

### Client testimonials

*"The Bioguard product, tested in its ready-to-use state, performed exceptionally well"*

**Law Laboratories results testing Disinfectant Cleaning Solution against MRSA**

*"Bioguard disinfectant shows bacterial activity to at least a dilution of 1:128 against the tested organism. As the working dilution of Bioguard is 1:15, this represents a total dilution of 1:1920 and is indicative of the efficiency of the product. I have no reservation about recommending Bioguard disinfectant for use in the American Mission Hospital"*

**Test conclusion, American Mission Hospital, Bahrain**

*"It's really good having a product you can use across a wide variety of surfaces, especially in my surgery where there's a lot of office furniture and fabrics"*

**Food Processor**

Ready to order or to ask a question?



PO Box 858, Northampton,  
NN2 1AF, United Kingdom

T: +44(0)1604 710333  
F: +44(0)1604 791587

info@bioguardhygiene.com  
www.bioguardhygiene.com

## Disinfectant cleaning solution

Non-hazardous, friendly to handle, easy on the environment and highly economic in use.

### Recommended environments and applications

Bioguard disinfectant cleaning solution is formulated for use in any healthcare environment where there's a risk of infection from dangerous micro-organisms.

These include:

- Hospitals
- Mortuaries
- First-aid posts
- Dental surgeries
- Food preparation areas
- GP surgeries
- Laboratories
- Emergency services
- Nursing & care homes
- Medical equipment, supplies & rental

Specific cleaning applications include:

- Work surfaces, floors and walls
- Fabric & Furnishings
- Equipment & Utensils
- Whole-room decontaminations



*"The deep cleaning exercise we carried out aboard ship in December, using your products, seems to have been very successful. The outbreak of viral gastroenteritis affecting the ship rapidly died down after she left Southampton"*

Multi-national cruise line operator

### Products Available

Dispensing option	Pack size	Product code
Trigger spray	12 x 500ml	BRTS1
1l Refill	6 X 1l	BR001
5l Refill	2 x 5l/4 x 5l	BR005

#### Concentrated solution

(Product dilution rates: heavy-duty cleaning 1:15; medium-duty cleaning 1:40; light-duty cleaning 1:75)

1l	12 x 1l	BC001
5l	2 x 5l/4 x 5l	BC005



Ready to order or to ask a question?



PO Box 858, Northampton,  
NN2 1AF, United Kingdom

T: +44(0)1604 710333  
F: +44(0)1604 791587

info@bioguardhygiene.com  
www.bioguardhygiene.com